


# The development of a robust methodology to evaluate inhalation capsule puncture performance

F. Díez<sup>2</sup>, B. Torrisi<sup>1</sup>, J.C. Birchall<sup>1</sup>, S.A. Coulman<sup>2</sup>, B. E. Jones<sup>1,2</sup>


<sup>1</sup> School of Pharmacy & Pharmaceutical Sciences, Cardiff University, Cardiff, CF10 3NB, UK, <sup>2</sup> Qualicaps Europe S.A.U., 28108 Alcobendas (Madrid), Spain

## PURPOSE

- Two-piece gelatin capsules have been used traditionally as an oral dosage form and as unit dose containers for a powdered drug for use in dry powder inhalers (DPI).
- Hypromellose is an alternative capsule material that has been shown to possess better functional properties than gelatin when used in DPI<sup>1,2</sup>.
- Both materials are currently used in DPI but a standardised methodology to evaluate capsule puncture performance, which is an essential property, does not exist.
- This study aimed to develop a robust methodology to determine potential differences in the puncture characteristics of different capsule materials, to assist in the development of hard capsules for this application.




**Figure 1.** Puncture pin from Plastiapipe Monodose Mod.7 inhaler




**Figure 2.** Zwitter Tester showing capsule holding device

**Figure 3.** Example of puncture holes in capsule caps conditioned at 33% RH: A, hypromellose; B, gelatin



## METHODS

- A steel conical tipped pin from a commercial DPI device (Plastiapipe S.p.a., Monodose Mod.7, 2 x 1 pin), see **Fig. 1**, was mounted in a bespoke miniaturised materials testing machine (Zwick® Testing Machines Ltd, UK), attached to an XForce P 500N load cell, see **Fig. 2**. The equipment is designed to measure small changes in force (accuracy  $\pm 1\%$  of the measured value) during a measurable displacement.
- A stainless steel bushing from a capsule-filling machine (Qualicaps), held a size 3 capsule in a fixed position directly below the steel pin.
- Hypromellose and gelatin capsules ( $n=10$  per each test), conditioned over saturated solutions of Calcium chloride (33% RH) and Lithium Chloride (11%RH) at 22°C for 1 week, to simulate poor storage conditions<sup>3</sup>.
- Capsules were punctured by the pin at a speed of 10 mm/s and the displacement of the pin and the resulting forces were recorded on a force-displacement curve<sup>4,5</sup>.
- Punctured capsules were subsequently removed from the bush for visual inspection, see **Fig. 3**.



**Figure 4.** Force(N)-Displacement(mm) curves for Hypromellose & Gelatin capsules conditioned at 11% and 33% RH

**Table 1.** Force/deformation values for Hypromellose and Gelatin capsules

Capsule	Hypromellose		Gelatin	
	11% RH	33% RH	11% RH	33% RH
Maximum force (N)	3.36 $\pm$ 0.35	3.43 $\pm$ 0.41	4.76 $\pm$ 0.62	4.45 $\pm$ 0.39
dL at F <sub>max</sub> (mm)	0.45 $\pm$ 0.07	0.49 $\pm$ 0.06	0.51 $\pm$ 0.03	0.51 $\pm$ 0.03

## RESULTS

- Repeated force-displacement profiles were highly reproducible for each of the capsule materials. However, HPMC and gelatin capsules possessed different signature profiles, characterised by differences in the penetration event, see **Figure 4**.
- Gelatin capsules: showed a rapid drop in force after puncture indicating that the pin had lost contact with the shell wall; the force then increased as the flap regained contact. This was more marked after lower RH conditioning indicating particles of the shell wall had become detached
- Hypromellose capsules: compared to gelatin the drop in force after puncture was less rapid and declined less; there was a much reduced second peak compared to gelatin. The difference probably being due to a lower elasticity of these puncture flaps compared to the gelatin ones.
- The force required to puncture hypromellose capsules was lower than for gelatin and occurred at a shorter deformation distance, see **Table 1**.

## CONCLUSIONS

- A rapid and robust methodology has been developed that is able to characterise penetration of a hard shell capsule by the pins that are employed in DPI.
- The sensitivity and reproducibility of the methodology enables users to describe differences in the capsule materials. This could have a significant impact on the design, development and quality assurance of hard shell capsules for use in DPI<sup>5</sup>.

## BIBLIOGRAPHY

- Jones, B.E., *Drug Deliv. Tech.*, 2003, **3**(6), 52-57, "Quali-V®-I: A new key for dry powder inhalers"
- Birchall, J.C., Jones, B.E., Morrissey, A., et al., *Drug Dev. Ind. Pharm.*, 2008, **34**, 870-876, "A comparison of the puncturing properties of gelatin and hypromellose capsules for use in dry powder inhalers"
- Renswouwe, D.C., van Laarhoven, A.C.M. et al., *J. Pharm. Pract.*, 2010, **23**, 548-552, "Storage conditions for inhalation capsules: consequences of incorrect storage and adherence in daily practice"
- Torrisi, B.M., Coulman, S. et al., Report on use of a Zwitter test device to measure inhalation capsule puncture, Cardiff University report prepared for Qualicaps Europe, March, 2013
- Torrisi, B.M., Birchall, J.C. et al., *Int. J. Pharm.*, 2013, **456**, 545-552, "The development of a sensitive methodology to characterise hard capsule shell puncture by dry powder inhaler pins"

We would like to thank José Luis Encinas of Qualicaps Europe, S.A.U. for supplying the capsule holder.